




## 2,5-Dimethoxyamphetamine

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.



### 1. GENERAL INFORMATION

<b>IUPAC Name:</b>	1-(2,5-dimethoxyphenyl)propan-2-amine
<b>CAS#:</b>	2801-68-5
<b>Synonyms:</b>	2,5-DMA; 2,5-dimethoxy- $\alpha$ -methyl-benzeneethanamine; 2,5-dimethoxy- $\alpha$ -methylphenethylamine
<b>Source:</b>	DEA Reference Material Collection
<b>Appearance:</b>	White powder (HCl and HBr)
<b>UV<sub>max</sub>(nm):</b>	Not Determined

### 2. CHEMICAL AND PHYSICAL DATA

#### 2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Melting Point (°C)
Base	C <sub>11</sub> H <sub>17</sub> NO <sub>2</sub>	195	Not Determined
HCl	C <sub>11</sub> H <sub>17</sub> NO <sub>2</sub> · HCl	231	108.8
HBr	C <sub>11</sub> H <sub>17</sub> NO <sub>2</sub> · HBr	276	121.5 and 130.4



## 2,5-Dimethoxyamphetamine

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.



### 3. QUALITATIVE DATA


#### 3.1 NUCLEAR MAGNETIC RESONANCE

##### Method NMR D<sub>2</sub>O

*Sample Preparation:* Dilute analyte to ~10 mg/mL in D<sub>2</sub>O containing TSP for 0 ppm reference and maleic acid as quantitative internal standard.

**Instrument:** 400 MHz NMR spectrometer  
**Parameters:** Spectral width: at least containing -3 ppm through 13 ppm  
Pulse angle: 90°  
Delay between pulses: 45 seconds

<sup>1</sup>H NMR: 2,5-dimethoxyamphetamine HCl Lot # AKB29A





## 2,5-Dimethoxyamphetamine

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.



### 3.2 Gas Chromatography/Mass Spectrometry

**Sample Preparation:** Dilute analyte ~1 mg/mL base extracted into chloroform.

**Instrument:** Agilent gas chromatograph operated in split mode with MS detector

**Column:** DB-1 MS (or equivalent); 30m x 0.25 mm x 0.25  $\mu$ m

**Carrier Gas:** Helium at 1 mL/min

**Temperatures:** Injector: 280°C

MSD transfer line: 280°C

MS Source: 230°C

MS Quad: 150°C

Oven program:

1) 100°C initial temperature for 1.0 min

2) Ramp to 300°C at 12 °C/min

3) Hold final temperature for 9.0 min

**Injection Parameters:** Split Ratio = 20:1, 1  $\mu$ L injected

**MS Parameters:** Mass scan range: 30-550 amu


Threshold: 100

Tune file: stune.u

Acquisition mode: scan

**Retention Time:** 8.077 min

EI Mass Spectrum: 2,5-dimethoxyamphetamine HCl Lot # AKB29A





# 2,5-Dimethoxyamphetamine

The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.





## 3.3 INFRARED SPECTROSCOPY (FTIR)

**Instrument:** FTIR with diamond ATR attachment (3 bounce)

**Scan Parameters:**  
Number of scans: 32  
Number of background scans: 32  
Resolution: 4 cm<sup>-1</sup>  
Sample gain: 8  
Aperture: 150

FTIR ATR (Diamond, 3 Bounce): 2,5-dimethoxyamphetamine HCl Lot # AKB29A



## **2,5-Dimethoxyamphetamine**

*The Drug Enforcement Administration's Special Testing and Research Laboratory generated this monograph using structurally confirmed reference material.*

### **4. ADDITIONAL RESOURCES**

[Forendex](#)

[Wikipedia](#)