

1. GENERAL INFORMATION

IUPAC Name:	3-(1,3-benzenodioxol-5-yl)-N,2-dimethylpropan-1-amine
CFR:	Not Scheduled (as of 2/2013)
CAS #:	N/A
Synonyms:	3,4-methylenedioxyamphetamine methylene homolog, 3,4-methylenedioxyamphetamine methyl homolog, Heliomethylamine
Source:	DEA Reference Material Collection
Appearance:	White powder (HCl)
Kovat's Index:	Pending
UV_{max} (nm):	Not Determined

2. CHEMICAL AND PHYSICAL DATA

2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Melting Point (°C)
Base	C ₁₂ H ₁₇ NO ₂	207	Not Determined
HCl	C ₁₂ H ₁₇ NO ₂ · HCl	243	144.9

3. ADDITIONAL RESOURCES

No resources identified as of 02/15/2013.

4. QUALITATIVE DATA

4.1 NUCLEAR MAGNETIC RESONANCE


Method NMR D₂O

Sample Preparation: Dilute analyte to ~5 mg/mL in D₂O containing TSP for 0 ppm reference and maleic acid as quantitative internal standard.


Instrument: Varian Mercury 400 MHz NMR spectrometer with proton detection probe

Parameters: Spectral width: at least containing -3 ppm through 13 ppm
Pulse angle: 90°
Delay between pulses: 45 seconds
Number of scans (NT): 8
Number of steady state scans: 0
Oversampling: 4 or more
Shimming: automatic gradient shimming of Z1-4 shims
Phasing, Drift Correction: automatic or manual


¹H NMR: 3,4-MDMA methylene homolog HCl; lot 0430048-14, D₂O, 400 MHz



^1H NMR: 3,4-MDMA methylene homolog HCl; lot 0430048-14, D_2O , 400 MHz



^1H NMR: 3,4-MDMA methylene homolog HCl; lot 0430048-14, 400 MHz



4.2 GAS CHROMATOGRAPHY/MASS SPECTROMETRY

Sample Preparation: Dilute analyte to ~4 mg/mL in chloroform.

Instrument: Agilent gas chromatograph operated in split mode with MS detector

Column: HP-1 MS; 30m x 0.25 mm x 0.25 μ m

Carrier Gas: Helium at 1 mL/min

Temperatures:
Injector: 280°C
MSD transfer line: 280°C
MS Source: 230°C
MS Quad: 150°C

Oven program:

1) 100°C initial temperature for 1.0 min

2) Ramp to 300°C at 12°C/min

3) Hold final temperature for 9.0 min

Injection Parameters: Split Ratio = 20:1, 1 μ L injected

MS Parameters: Mass scan range: 30-550 amu


Threshold: 90

Tune file: stune.u

Acquisition mode: scan

Retention Time: 9.471 min


GC-MS: 3,4-MDMA methylene homolog HCl; lot 0430048-14



4.3 INFRARED SPECTROSCOPY (FTIR)

Instrument: FTIR with diamond ATR attachment (3 bounce)
Scan Parameters: Number of scans: 32
Number of background scans: 32
Resolution: 4cm^{-1}
Sample gain: 8
Aperture: 150

FTIR ATR (Diamond, 3 bounce): 3,4-MDMA methylene homolog HCl; lot 0430048-14



FTIR ATR (Diamond, 3 bounce): 3,4-MDMA methylene homolog HCl; lot 0430048-14

