

Analytical Profile of 30C-NBOMe

Michael Martin, Danielle St. Germaine, Cayman Chemical Company, Ann Arbor, MI 48108

Erin Tracy, Patrick Chepaitis, Georgia Bureau of Investigation, Headquarters Laboratory –
Decatur, GA

Table of contents:

Part 1: Cayman Chemical Data:

[Figure 1.](#) 30C-NBOMe GC/MS Data

[Figure 2.](#) 30C-NBOMe 1H NMR Data

[Figure 3.](#) 30C-NBOMe 1H NMR Data, Enlarged

[Figure 4.](#) 30C-NBOMe ESI-MS Data

Part 2: Georgia Bureau of Investigation Data:

[Figure 5.](#) 30C-NBOMe GC/MS Data

[Figure 6.](#) 30C-NBOMe Expanded GC/MS Data

[Figure 7.](#) 30C-NBOMe FTIR Data

[Figure 8.](#) 30C-NBOMe Expanded FTIR Data

[Figure 9.](#) 30C-NBOMe HPLC Data

[Figure 10.](#) 30C-NBOMe UV Data

Part 3: External Link

Part 1. Cayman Chemical Company Data

Name: 30C-NBOMe

Synonyms: 2-(4-chloro-2,5-dimethoxyphenyl)-N-(3,4,5-trimethoxybenzyl)ethanamine

CAS#: N/A


Molecular Formula: C₂₀H₂₆ClNO₅

Molecular Weight: 395.5 g/mol

SMILES: COc2cc(CCNc1cc(OC)c(OC)c(OC)c1)c(cc2Cl)OC

InChI Key: ZQYYVTABADQBTJ-UHFFFAOYAT

InChI: InChI=1/C20H26ClNO5/c1-23-16-11-15(21)17(24-2)10-14(16)6-7-22-12-13-8-18(25-3)20(27-5)19(9-13)26-4/h8-11,22H,6-7,12H2,1-5H3



Background:

30C-NBOMe is a N-alkylated analogue of the phenethylamine 2C-C, a known hallucinogen that stimulates monoamine receptor activity and inhibits the reuptake of serotonin (IC₅₀ = 31uM) and norepinephrine (IC₅₀ = 63uM) in rat brain synaptosomes.^{1,2} 30C-NBOMe exhibits a 3,4,5 tri-methoxy benzyl group bridged by the nitrogen to the 2C-C core. The aryl benzyl moiety has been shown to increase affinity and selectivity for the 5-HT_{2A} receptor in the non-chlorinated 2C-H phenethylamine series³. Moreover, the 3,4,5 tri-methoxy benzyl group resembles the hallucinogenic alkaloid mescaline, less one carbon in the alkyl chain. The physiological and toxicological properties of this compound are not known.

1. Nagai, F., Nonaka, R., Satoh Hisashi Kamimura, K. *Eur. J. Pharmacol.* **2007**. 559(2-3), 132-137.
2. Nonaka, R., Nagai, R., Ogata, A., *et al. Biol. Pharm. Bull.* **2007**. 30(12), 2328-2333.
3. Braden, M.R., Parrish, J.C., Naylor, J.C., Nichols, D.E. *Mol Pharm.* **2006**. 70(6), 1956-1964.

Gas chromatography/Mass spectrometry:

Experiment Parameters: **Instrument:** Agilent 6890 GC / 5973 MSD

Column: 30mx0.32mm, 0.5um Rtx-5MS

Carrier Gas: Helium Flow: 2mL/min

Inlet temp: 300 °C, 15:1 split

Oven Program: Initial temp: 240 °C, Ramp to 300 °C at 30C/min,
Hold at 300 °C for 27 minutes

Transfer Line Temp: 300 °C

MS Source: 230°C

MS Quad: 150°C

Mass Scan Range: 40-600 amu

Threshold: 150

Tune File: stune.u

Figure 1: Gas Chromatography/Mass Spectrometry of 3OC-NBOMe (hydrochloride)




Figure 2. ¹H Nuclear Magnetic Resonance Spectroscopy of 3OC-NBOMe (hydrochloride)




Figure 3. ¹H Nuclear Magnetic Resonance Spectroscopy of 30C-NBOMe (hydrochloride), Enhanced for Detail





Figure 4. ESI-MS of 30C-NBOMe (hydrochloride)

Experiment Parameters:

Instrument: Thermo LCQ Ion Trap MS
 Sample Infusion: 500 µl/min
 Mode: +ESI
 Mass range: 100-700 amu
 Capillary Temperature: 275 C



Part 2. Georgia Bureau of Investigation Data

Name: 30C-NBOMe

Synonyms: 2-(4-chloro-2,5-dimethoxyphenyl)-N-(3,4,5-trimethoxybenzyl)ethanamine

CAS#: N/A

Molecular Formula: C₂₀H₂₆ClNO₅

Molecular Weight: 395.5 g/mol


Source: Cayman Chemicals lot # 0455137-3

Appearance: white, crystalline solid

SMILES: COc2cc(CCNCc1cc(OC)c(OC)c(OC)c1)c(cc2Cl)OC

InChI Key: ZQYYVTABADQBTJ-UHFFFAOYAT

InChI: InChI=1/C20H26ClNO5/c1-23-16-11-15(21)17(24-2)10-14(16)6-7-22-12-13-8-18(25-3)20(27-5)19(9-13)26-4/h8-11,22H,6-7,12H2,1-5H3



Gas Chromatography/Mass Spectrometry (GCMS):

Sample Preparation: ~ 2mg/mL in EtOH

Instrument: Agilent 7890A gas chromatograph / 5975C mass spectrometer

GC Parameters **Column:** HP Ultra-1 12m × 200µm × 0.33µm

Carrier gas: Helium

Oven program: 1) initial temperature 60°C for 2 min.

 2) ramp to 275°C at 35°C/min.

 3) hold at 275°C for 12 min.

Injection parameters: 1µL injection; pulsed splitless

Inlet temperature: 270°C

Constant Pressure: 8 psi

MS Parameters **Transfer Line temperature:** 280°C

Source temperature: 230°C

Quad temperature: 150°C

Mass scan range: 40-500 amu

Tune file: stune

Figure 5: Gas Chromatography/Mass Spectrometry of 30C-NBOMe (hydrochloride)





Figure 6: Expanded Gas Chromatography/Mass Spectrometry of 3OC-NBOMe (hydrochloride)



Infrared Spectroscopy (FTIR):

Instrument: Nicolet 4700 with diamond ATR attachment (3 bounce)

Scan Parameters: Number of scans: 32

Number of background scans: 32

Resolution: 4 cm^{-1}

Sample gain: 8

Figure 7: Infrared Spectroscopy of 30C-NBOMe (hydrochloride)





Figure 8: Expanded Infrared Spectroscopy of 3OC-NBOMe (hydrochloride)




High Performance Liquid Chromatography (HPLC):**Instrument:** Agilent 1260 Infinity**LC Parameters:** **Column:** Zorbax Eclipse Plus C18 3.0 x 100mm 35-micron**Column compartment:** 40°C**Mobile phase:**A: Ammonium Phosphate buffer

B: Acetonitrile

<u>Time</u>	<u>%A</u>	<u>%B</u>	<u>Flow</u>
0	95	5	1.000 mL/min
0.67	95	5	1.000 mL/min
0.68	85	15	1.000 mL/min
2.00	85	15	1.000 mL/min
6.70	15	85	1.000 mL/min
6.70	15	85	1.000 mL/min
7.50	15	85	1.000 mL/min
8.00	5	95	1.000 mL/min
8.00	5	95	2.700 mL/min
8.50	5	95	2.700 mL/min
9.75	95	5	1.000 mL/min
10.25	95	5	1.000 mL/min

DAD Parameters: **Scan:** 220 – 340nm**Step:** 1.0nm**Slitwidth:** 4nm

Figure 9: High Performance Liquid Chromatography of 30C-NBOMe (hydrochloride)



Ultraviolet Spectrophotometry (UV):

Instrument: HP 8453

Parameters: **Source:** deuterium lamp

Path length: 1cm


Transmittance

Wavelength range: 220 – 340nm

Integration time: 0.5 sec

Interval: 1nm

Figure 10: Ultraviolet Spectrophotometry of 30C-NBOMe (hydrochloride)



Part 3. External Links:

[Forendex](#)