

## 1. GENERAL INFORMATION

<b>IUPAC Name:</b>	naphthalen-1-yl(1-pentyl-1 <i>H</i> -indol-3-yl)methanone
<b>CFR:</b>	Schedule I
<b>CAS #:</b>	209414-07-3
<b>Synonyms:</b>	AM678, 1-pentyl-3-(1-naphthoyl)indole
<b>Source:</b>	DEA Reference Material Collection
<b>Appearance:</b>	Off-white powder
<b>Kovat's Index:</b>	Pending
<b>UV<sub>max</sub>:</b>	218.1, 314.6 nm

## 2. CHEMICAL AND PHYSICAL DATA

### 2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Melting Point (°C)
Base	C <sub>24</sub> H <sub>23</sub> NO	341	51.9

### 3. ADDITIONAL RESOURCES

Kneisel S, Westphal F, Rosner P, et.al. Cannabinoidmimetika: Massenspektren und IR-ATR-Spektren neuer Verbindungen aus den Jahren 2009/2010. *Toxichem Krimtech.* 2011; 78(1):23-35.

[Forendex](#)

[Wikipedia](#)

### 4. QUALITATIVE DATA

#### 4.1 NUCLEAR MAGNETIC RESONANCE

##### *Method NMR CDCl<sub>3</sub>*


*Sample Preparation:* Dilute analyte to ~20 mg/mL in deuteriochloroform (CDCl<sub>3</sub>) containing TMS for 0 ppm reference and dimethylfumarate as quantitative internal standard.

***Instrument:*** Varian Mercury 400 MHz NMR spectrometer with proton detection probe


***Parameters:***

- Spectral width: at least containing -3 ppm through 13 ppm
- Pulse angle: 90°
- Delay between pulses: 45 seconds
- Number of scans (NT): 8
- Number of steady state scans: 0
- Oversampling: 4 or more
- Shimming: automatic gradient shimming of Z1-4 shims
- Phasing, Drift Correction: automatic or manual


1H NMR: JWH-018 Lot # ALB045RC/183-1; CDCl<sub>3</sub>; 400MHz




1H NMR: JWH-018 Lot # ALB045RC/183-1; CDCl<sub>3</sub>; 400MHz



1H NMR: JWH-018 Lot # ALB045RC/183-1; CDCl<sub>3</sub>; 400MHz



1H NMR: JWH-018 Lot # ALB045RC/183-1; CDCl<sub>3</sub>; 400MHz



## 4.2 GAS CHROMATOGRAPHY/MASS SPECTROMETRY

*Sample Preparation:* Dilute analyte to ~1 mg/mL in CHCl<sub>3</sub>.

**Instrument:** Agilent gas chromatograph operated in split mode with MS detector

**Column:** DB-1 MS or equivalent; 30m x .25mm x .25µm

**Carrier Gas:** Helium at 1 mL/min

**Temperatures:** Injector: 280°C

MSD transfer line: 280°C

MS Source: 230°C

MS Quad: 150°C

Oven program:

1) 90°C initial temperature for 2.0 min

2) Ramp to 300°C at 14°C/min

3) Hold final temperature for 10.0 min

**Injection Parameters:** Split Ratio = 25:1, 1 µL injected

**MS Parameters:** Mass scan range: 34-550 amu


Threshold: 100

Tune file: stune.u

Acquisition mode: scan

**Retention Time:** 20.011 minutes


EI Mass Spectrum: JWH-018 Lot # ALB045RC/183-1



### 4.3 INFRARED SPECTROSCOPY (FTIR)

**Instrument:** FTIR with diamond ATR attachment (3 bounce)  
**Scan Parameters:** Number of scans: 32  
Number of background scans: 32  
Resolution:  $4\text{cm}^{-1}$   
Sample gain: 8  
Aperture: 150

FTIR ATR (Diamond, 3 Bounce): JWH-018 Lot # ALB045RC/183-1



FTIR ATR (Diamond, 3 Bounce): JWH-018 Lot # ALB045RC/183-1

